lumie[®]

bodyclock[™]

luxe 750D

Safety

If the product is damaged please contact us before use.

Keep the unit away from water and damp. 

For indoor use only. 

The lamp cover gets warm when the light is on, so keep the unit clear of anything that may be damaged by heat.

Use only the mains power adaptor supplied.

Keep the power supply and cable away from heated surfaces. If either are damaged they should be replaced with a Lumie approved product.

This appliance does not have an on/off switch. To power down or reset the appliance, remove the mains power adaptor from the wall socket.

Do not attempt to service any part of the unit.

Not to be used by or on individuals who are unconscious, anaesthetised or incapacitated e.g. under the influence of alcohol.

Electromagnetic compatibility

This product conforms to the electromagnetic compatibility (EMC) standard for medical devices and is unlikely to cause interference affecting other electrical and electronic products in its vicinity. However, radio transmitting equipment, mobile phones and other electrical/electronic equipment can interfere with this product; if this happens simply move the product away from the source of interference. In rare cases where product function is lost due to EMC effects, disconnect the product from the power source for 5 seconds and then restart. If the problem persists, please call the Lumie Careline.

For further detailed information on the electromagnetic environment in which the device can be operated refer to the website: www.lumie.com/pages/e-m-c.

Disposal

At the end of its life, do not throw the unit out with normal domestic waste.  Please take it to an official recycling point to be disposed of responsibly.

 SN	Serial number		Do not place in household waste. Dispose of at a recycling centre.		Double-insulated (Class II electrical)
	European Declaration of Conformity: Class I Medical Device		Read instructions for use		AC (alternating current) power supply
	Manufacturer: Lumie, 3 The Links, Bar Hill, Cambridge, CB23 8UD, UK.		Caution		Member of the 'Green Dot' recycling scheme
	Keep dry		For indoor use only		Recyclable

Introduction

Bodyclock Luxe is a wake-up light and music player with extra bedtime features to help you fall asleep.

Wake-up light

Bodyclock Luxe wakes you with a sunrise that naturally helps to regulate your sleep and wake patterns. You'll wake up in a better mood, find it easier to get up and feel more energetic.

Ideally, you should use the wake-up light all year round to help keep your sleep and wake cycle on track.

In clinical trials, 'dawn simulators' such as Bodyclock Luxe have been shown to improve mood, energy, productivity and the quality of sleep and awakening¹. They can also help sufferers of SAD (Seasonal Affective Disorder) and winter blues cope with dark winter mornings².

Bedside lighting

Bodyclock Luxe provides adjustable bedside lighting and a relaxing sunset to help you unwind at bedtime. Use together with the low-blue light feature to minimise the alerting effect of blue light and prepare for sleep³.

Music player

Bodyclock Luxe has DAB radio and can also play the radio and music wirelessly from a smartphone or computer via Bluetooth, from a USB flash drive or from a music player via an audio cable.

1. Acta Psychiatr Scand. 1993 Jul;88(1):67-71.

2. J Affect Disord. 2015 Jul 15;180:87-9.

3. J Pineal Res. 2012 Aug;53(1):47-59.


Getting started

Bodyclock Luxe is supplied with the following options already set to get you started:

- A daily alarm at 07:00.
- A 30-minute sunrise reaching 80% brightness to wake you naturally.
- A 30-minute sunset with gradually fading light to help you fall asleep.
- A display that's dim during the night and brighter during the day.

Here's how to set up Bodyclock Luxe for the first time.

To set the time

- First, connect the mains power adaptor to the socket on the base of the unit, and feed the power cable through the channel provided.
- Plug in the mains power adaptor and switch on at the mains.

A  cursor will prompt you to enter the essential settings:

Turn the right dial to choose a language.

Press to confirm.

The cursor moves to the next setting.

Turn the dial to set the day.

Press to confirm.

Now edit the time.

Turn the dial to set the hour and press to confirm.


Do the same for the minutes.

Bodyclock Luxe is now ready to use as a wake-up light, bedside light and bedside music player.

Using Bodyclock Luxe

Front panel

You operate Bodyclock Luxe using the buttons and dials on the front panel.

The four buttons perform specific functions indicated by the icon next to them on the display. The bottom right icon indicates which mode you're in.

In **Light Mode** the bottom right shows a greyed-out audio icon or 

In this mode the two dials adjust the brightness of the light and the display.

In **Audio Mode** the bottom right shows an audio icon:


In this mode the two dials adjust the volume and allow you to select audio tracks.

The five preset buttons allow you to save and play your favourite radio stations.

Back panel

The back panel allows you to connect other devices to Bodyclock Luxe:

Using the alarm

Bodyclock Luxe has two types of alarm:

DAILY gives you a sunrise alarm at the same time every day. You choose each evening whether you want to activate the alarm for the following morning.

WEEKLY allows you to set the alarm time for each day individually. For example, you may want to wake up earlier on Monday and have no alarm at the weekend.

To choose DAILY or WEEKLY alarm type

- Press the settings button .

The display shows the first settings menu, ALARM:

- Press the right dial to edit ALARM TYPE.

The  cursor changes to  to show that you are editing the item.

- Turn to choose DAILY or WEEKLY and press to confirm.

The cursor will move to the next setting.

To set up a daily alarm

If you selected DAILY you need to set an alarm time:

- Press the right dial to edit the hour.
- Turn to set the hour and press to confirm.
- Do the same for the minutes.
- Press the button next to  to exit from settings.

After setting, the next alarm is automatically activated and is shown on the top line of the display:

To set up weekly alarms

If you selected WEEKLY you can specify for each day whether the alarm is OFF or ON, and if ON, the wake-up time:

- Press the right dial to edit the day; for example, MON.

The  cursor changes to  to show that you are editing the item.

- Turn to choose OFF or ON and press to confirm.

If you selected ON, an additional SET TIME item appears below the day.

- Turn to set the hour and press to confirm.
- Do the same for the minutes.
- Repeat for each day.
- Press the button next to  to exit from settings.

After setting, the next alarm is automatically activated and is shown on the top line of the display:

You are now ready to use Bodyclock Luxe with the default sunset and sunrise settings.

To change these settings or choose a wake-up sound see *Choosing sunset and sunrise options*, page 18.

Using the bedside light

Bodyclock Luxe can be used as a bedside light that you can set to your preferred brightness.

In **Light Mode** use the two dials to adjust the brightness of the light and display.

In **Light Mode** the bottom right audio icon is greyed-out or 

If you're in **Audio Mode** press the light button next to  to switch to **Light Mode**.

To use the bedside light

- Turn the left dial to adjust the brightness of the light, shown from 00 to 20 on the left of the display.

To switch to low-blue light

The blue part of the spectrum is known to have an alerting effect so it's good to minimise this at bedtime. Choose low-blue light (this looks warmer and more yellow than the white light) to optimise your bedroom lighting for relaxing sleep.

- Press the left dial to switch between normal white light and low-blue light.

To adjust the display brightness

- Turn the right dial to adjust the brightness of the display, shown from 00 to 20 on the right of the display.

If you have dimmed the display, and want to turn it on at night to see the time, press the right dial.

Note: The display brightness can't be manually adjusted during a sunset.


Using the audio

Bodyclock Luxe incorporates a high-quality audio system with DAB radio that can also play music from a smartphone or other device via Bluetooth, from a USB flash drive or from the auxiliary input.

You can listen through the built-in stereo speakers or plug in earphones or headphones for private listening.

In **Audio Mode** use the two dials to adjust the volume and select tracks.

In **Audio Mode** the bottom right shows an audio icon. If you're in **Light Mode** press the audio button to switch to **Audio Mode**.

To select your audio source

- If audio is currently playing, first turn it off by pressing the audio button.
- Press the audio button next to  to step through the audio options:

DAB radio	Bluetooth	USB	Auxiliary input	Internal sounds
DAB			AUX	

After two seconds the audio source shown is selected.

- Adjust the volume with the left dial.

The following sections give more information about each of the audio options.

To use DAB radio

DAB provides high-quality stereo sound and a wide range of stations.

- Press the audio button until the icon next to it shows **DAB**.
- Turn the right dial to select the station.

The station is shown on the bottom line of the display:

The signal strength is shown by the  icon. If necessary adjust the telescopic aerial for best reception.

To save a station as a preset

You can save up to five DAB for instant listening:

- Turn the right dial to select the station you want to save.
- Press and hold a preset button until the display shows **SAVED**.
- Press the preset key to listen to your saved station.

To play the radio or music via Bluetooth

You can listen to the radio or play your own music using Bodyclock Luxe connected to your smartphone, tablet or other device via Bluetooth.

- Press the audio button until the icon next to it shows .
- On your other device, search for available Bluetooth devices and select **Lumie 700** for pairing.

The display will show **BT PAIRED** and you can now play music through Bodyclock Luxe:

To play music from a USB flash drive

You can play music from a library of audio files on a USB flash drive. Bodyclock Luxe supports MP3 files using one level of folders.

Note: The USB socket should only be connected to a USB flash drive, not an external hard drive. It should not contain more than 200 folders or 2,000 music files. The drive should be formatted as FAT32 and have a maximum size of 32 GB. Lumie recommends Sandisk USB drives.

- Connect the USB drive to the USB socket on the back panel.
- Press the audio button until the icon next to it shows .

The display shows the first folder on the USB drive and a list of files. If you have a large number of music files there may be a delay before the list appears.

- Turn the left dial to step between folders on the USB drive.
- Turn the right dial to highlight the track you want to listen to and press to select it.

Use the right dial to skip between tracks during play.

To play music from the AUX input

You can play music from any device with an audio output, such as a smartphone or MP3 player.

- Connect the device to the AUX socket on the back panel using a cable terminated in a mini-jack plug.

IMPORTANT: Only connect the AUX socket to the headphone socket of a smartphone or MP3 player.

- Press the audio button until the icon next to it shows **AUX**:

You can either adjust the sound level on your music player or using the left dial on Bodyclock Luxe.


To turn off the audio

- Press the audio button so the icon next to it shows .

Using the sunset and sunrise

You can use Bodyclock Luxe to go to sleep with a sunset and wake up to a sunrise.

At bedtime

If you've selected DAILY you need to activate the sunrise alarm for the following morning:

- Press the alarm button 

The top line of the display shows the alarm time; for example:

If you've selected WEEKLY the alarm is activated automatically according to the schedule you've programmed.

To activate the sunset

When you're ready to go to sleep you can activate the sunset. Start the sunset with the light at full brightness if you like to read for a while first or low for a gently fading background light.

- Make sure you're in **Light Mode**.
- Turn the left dial to adjust the light. Press the light button  to start the sunset.

The light will gradually fade to help you fall asleep according to the sunset settings you have selected. You can choose to have the audio option you selected turn off or stay on at the end of the sunset, or gradually fade to off with the sunset; see *To set the sunset options*, page 18.

In the morning

You'll be woken up with a gradual sunrise according to the sunrise settings you have selected. You can choose to have a wake-up sound, play the radio or play a USB track at the alarm time; see *To set the sunrise options*, page 19.

To snooze

- Press any preset key.

This gives you another nine minutes before the sunrise (and audible alarm if you've selected one) is repeated.

To turn off the sunrise alarm

- Press the alarm button .


This also turns off the audible alarm if you've selected one.

- Turn the left dial to turn off the light.

Choosing sunset and sunrise options

To set the sunset options

- Press the settings button  to display the settings menus.
- Turn the left dial until the display shows the SUNSET menu:

- Press the right dial to edit the DURATION item.
- Turn to set the duration of the sunset and press to confirm.

Note: The actual duration may be less than this if you start the sunset with the light already dimmed.

- Press to edit the NIGHTLIGHT item.

This is the brightness of the light during the night. As you edit it the light will change to show the brightness you are setting.

- Turn to set a level from 1 to 10 or OFF and press to confirm.
- Press to edit the DISPLAY item.

This is the brightness of the display during the night. As you edit it the display will change to show the brightness you are setting.

- Turn to set a level from 1 to 20 or OFF and press to confirm.
- Press to edit the SOUND item.
- Turn to choose between OFF, MANUAL or SOUND and press to confirm.

OFF turns off the audio when the sunset starts, MANUAL leaves your current audio selection playing when you start the sunset, and SOUND allows you to select a soothing sleep sound to play during the sunset.

If you selected SOUND an additional TRACK item allows you to select the sleep sound you want to play.


If you selected MANUAL or SOUND an additional VOLUME item allows

you to specify how the volume should change during the sunset. You can choose OFF or ON, to have the sound turn off or stay on at the end of the sunset, or FADE, to have the sound fade to off with the sunset.

- Press the button next to  to exit from settings.

To set the sunrise options

- Press the settings button  to display the settings menus.
- Turn the left dial until the display shows the SUNRISE menu:

- Press the right dial to edit the DURATION item.
- Turn to set the duration of the sunrise and press to confirm.

If you do not want a sunrise set the DURATION to OFF. The LIGHT LEVEL and LOW BLUE items will then disappear from the menu.

- Press to edit the LIGHT LEVEL item.

This is the final brightness of the sunrise. As you edit it the light will change to show the brightness you are setting.

- Turn to set a level from 1 to 20 and press to confirm.
- Press to edit the LOW BLUE item.
- Turn to choose ON or OFF and press to confirm.

The default setting is OFF which gives you a bright sunrise to help you wake up. Choose ON if you want the low-blue sunrise which is warmer and more yellow.

- Press to edit the DISPLAY item.

This is the brightness of the display at the end of the sunrise. As you edit it the display will change to show the brightness you are setting.

- Turn to set a level from 1 to 20 or OFF to leave the display off, and press to confirm.

Note: Settings 0 to 3 will be too dim to view in daylight.

- Press to edit the ALARM SRC (alarm source) item.
- Turn to choose OFF, RADIO, SOUND or USB and press to confirm it.

The default setting is OFF, no audible alarm. Alternatively you can choose RADIO to wake up to a radio station, SOUND to play one of the internal wake-up sounds or USB to play an audio file in a folder on a USB flash drive.

If you selected SOUND or USB an additional TRACK item allows you to select the wake-up sound you want to play.

If you selected RADIO, SOUND or USB an additional VOLUME item allows you to specify the volume during the sunrise from 1 to 20. The volume will change accordingly as you change the setting.

- Press the button next to  to exit from settings.

Customisation options

This section summarises the alarm, sunrise, sunset and other options available with Bodyclock Luxe.

To access the settings menus

- Press the settings button .

The display shows the first settings menu, ALARM:

- Turn the left dial to step through the different menus.
- Turn the right dial to browse the current menu.

The currently highlighted item is shown by a  cursor to its right.

Longer menus will scroll to show you all the items.

To change a setting

- Press the right dial.

The  cursor changes to  to show that you are editing the item.

- Turn the right dial to change the item.
- Press the right dial to confirm the change.

The  cursor will move to the next item on the menu.

To exit from the settings menus

- Press the button next to .

The following table summarises each menu and menu item:

Menu	Menu items	Options	Default
ALARM	ALARM TYPE	DAILY, WEEKLY	DAILY
	MON to SUN	OFF, ON	ON
	SET TIME	HH:MM	07:00
SUNRISE	DURATION	OFF or 15 to 90 MIN in 5 minute steps	30 MIN
	LIGHT LEVEL	1 to 20	16
	LOW BLUE	OFF, ON	OFF
	DISPLAY	OFF or 1 to 20	16
	ALARM SRC	Alarm source: OFF, RADIO, SOUND or USB	OFF
	STATION *	P1 to P5	
	TRACK *	Filename	
VOLUME *	1 to 20	10	
SUNSET	DURATION	15 to 90 MIN in 5 minute steps	30 MIN
	NIGHTLIGHT	OFF or 1 to 10	OFF
	DISPLAY	OFF or 1 to 20	3
	SOUND	OFF, MANUAL or SOUND	MANUAL
	TRACK *	Filename	
	VOLUME *	OFF, ON or FADE	OFF
GENERAL	LANGUAGE	ENG, FRA, ESP, DEU, NED, SVE, ITA, NOR or DAN	ENG
	DAY	MON to SUN	MON
	TIME		13:00
DAB	AUTOSCAN	START	

* Displayed depending on previous settings.

Frequently asked questions

I tried out the sunrise but the light comes up really quickly – what's wrong?

You probably didn't allow enough time for a full sunrise. For example, the default sunrise takes 30 minutes so if the clock shows the time as 11:00 then you need to set the wake-up time to at least 11:30.

What happens when I unplug my Bodyclock Luxe or there's a power cut?

Your Bodyclock Luxe keeps time for 30 minutes; other settings are permanent, so you won't need to reprogram the unit, just plug it back in. If you need to adjust the time, see *TIME* in the *GENERAL* menu, page 22.

The display panel is too bright when I'm trying to sleep; how can I change it?

See *DISPLAY* in the *SUNSET* menu, page 22.

How do I restore the default factory settings?

Turn off the power at the mains. Press and hold preset key 3 as you turn the power back on. Release when the software version number is displayed.

How can I control the playing order of files on a USB flash drive?

The files are played in alphabetical order, so you can prefix each filename with a number such as 01, 02, 03 etc to specify the playing order. The drive should not contain more than 2,000 music files or 200 folders.

Why can't I find any DAB stations?

You may need to perform a scan to find the stations available in your current area.

Highlight *AUTO SCAN* in the *DAB* menu, and select *START*; see *AUTO SCAN* in the *DAB* menu, page 22.

Maintenance

Clean only with a dry or lightly moistened cloth.

Replacement mains power adaptors are available from Lumie.

Warranty

This unit is covered under warranty against any manufacturing defect for three years from the date of purchase. This is in addition to your statutory rights. If possible, use the original packaging (including internal pieces) to return your Lumie for service. Make sure the mains power adaptor is prevented from moving in transit and damaging the unit. If the unit is received damaged, through misuse or accidental damage, then we may have to charge to repair it; we will contact you first. Please call Lumie for further advice.

Technical description

Subject to change without notice.

Lumie Bodyclock Luxe 750D

100-240 V~, 50-60 Hz Lumie mains power adaptor (supplied) 

Operating temperature: 5 °C to 35 °C

Storage/transport temperature range: -25 °C to +70 °C

Class I Medical Device 

Contact us

Lumie is Europe's leading specialist in light therapy. Our products are based on the latest scientific studies and we focus solely on light therapy and its benefits.

If you have any suggestions that you feel could help us improve our products or service, or if there are problems you think we should know about, please call the Lumie Careline on +44 (0)1954 780500 or email info@lumie.com.

Lumie, 3 The Links, Trafalgar Way, Bar Hill, Cambridge, CB23 8UD, UK.

Lumie is a registered trademark of:

Outside In (Cambridge) Limited, incorporated in England and Wales.

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© 2017
Outside In (Cambridge) Limited, eingetragen im England und Wales.
Outside In (Cambridge) Ltd